

# ExpressionAtlas package vignette

*Maria Keays*

*2016-01-20*

## Expression Atlas

The [EMBL-EBI Expression Atlas](#) consists of hand-picked high quality datasets from [ArrayExpress](#) that have been manually curated and re-analyzed via the Expression Atlas analysis pipeline. The Expression Atlas website allows users to search these datasets for genes and/or experimental conditions, to discover which genes are expressed in which tissues, cell types, developmental stages, and hundreds of other experimental conditions.

The *ExpressionAtlas* R package enables download of pre-packaged data from Expression Atlas directly into an R session. Raw counts are provided for RNA-seq datasets, while normalized intensities are available for microarray experiments. Protocols describing how the data was generated are contained within the downloaded R objects, with more detailed information available on the [Expression Atlas website](#). Sample annotations are also included in the R object.

## Searching and downloading Expression Atlas data

### Searching

You can search for experiments in Atlas using the `searchAtlasExperiments()` function. This function returns a *DataFrame* (see [S4Vectors](#)) containing the results of your search. The first argument to `searchAtlasExperiments()` should be a character vector of sample properties, e.g. biological sample attributes and/or experimental treatments. You may also optionally provide a species to limit your search to, as a second argument.

```
suppressMessages( library( ExpressionAtlas ) )
atlasRes <- searchAtlasExperiments( "salt", "Oryza sativa" )
```

```
## Searching for Expression Atlas experiments matching your query
## http://www.ebi.ac.uk/arrayexpress/xml/v2/experiments?keywords=salt&gxa=TRUE&species=Oryza%20sativa .
## Query successful.
## Found 3 experiments matching your query.
```

```
atlasRes
```

```
## DataFrame with 3 rows and 4 columns
##      Accession                      Species
##      <character>                    <character>
## 1 E-GEOD-11175 Oryza sativa Japonica group
## 2 E-MTAB-1625  Oryza sativa Japonica Group
## 3 E-MTAB-1624  Oryza sativa Japonica Group
##                                     Type
##                                     <character>
## 1 transcription profiling by array
## 2 RNA-seq of coding RNA
```

```
## 3 transcription profiling by array
##
##
## 1           Transcription profiling by array of wild type and drought and salt tolerance
## 2           RNA-seq of coding RNA in rice (Oryza sativa L. Japonica nipponbare) under salt stress o
## 3 Transcription profiling by array of rice (Oryza sativa L. Japonica Nipponbare) under salt stress o
```

The *Accession* column contains the ArrayExpress accession of each dataset – the unique identifier assigned to it. The species, experiment type (e.g. microarray or RNA-seq), and title of each dataset are also listed.

## Downloading the data

To download the data for any/all of the experiments in your results, you can use the function `getAtlasData()`. This function accepts a vector of ArrayExpress accessions. The data is downloaded into a *SimpleList* object (see package [S4Vectors](#)), with one entry per experiment, listed by accession.

For example, to download all the datasets in your results:

```
allExps <- getAtlasData( atlasRes$Accession )
```

```
## Downloading Expression Atlas experiment summary from:
## http://wwwdev.ebi.ac.uk/gxa/experiments/E-GEOD-11175/E-GEOD-11175-atlasExperimentSummary.Rdata ...
## Successfully downloaded experiment summary object for E-GEOD-11175
## Downloading Expression Atlas experiment summary from:
## http://wwwdev.ebi.ac.uk/gxa/experiments/E-MTAB-1625/E-MTAB-1625-atlasExperimentSummary.Rdata ...
## Successfully downloaded experiment summary object for E-MTAB-1625
## Downloading Expression Atlas experiment summary from:
## http://wwwdev.ebi.ac.uk/gxa/experiments/E-MTAB-1624/E-MTAB-1624-atlasExperimentSummary.Rdata ...
## Successfully downloaded experiment summary object for E-MTAB-1624
```

```
allExps
```

```
## List of length 3
## names(3): E-GEOD-11175 E-MTAB-1625 E-MTAB-1624
```

To only download the RNA-seq experiment(s):

```
rnaseqExps <- getAtlasData(
  atlasRes$Accession[
    grep(
      "rna-seq",
      atlasRes$Type,
      ignore.case = TRUE
    )
  ]
)
```

```
## Downloading Expression Atlas experiment summary from:
## http://wwwdev.ebi.ac.uk/gxa/experiments/E-MTAB-1625/E-MTAB-1625-atlasExperimentSummary.Rdata ...
## Successfully downloaded experiment summary object for E-MTAB-1625
```

```
rnaseqExps
```

```
## List of length 1  
## names(1): E-MTAB-1625
```

To access an experiment summary, use the accession:

```
mtab1624 <- allExps[[ "E-MTAB-1624" ]]  
mtab1625 <- allExps[[ "E-MTAB-1625" ]]
```

Each dataset is also represented by a *SimpleList*, with one entry per platform used in the experiment. For RNA-seq data there will only ever be one entry, named **rnaseq**. For microarray data, there is one entry per array design used, listed by ArrayExpress array design accession (see below).

### RNA-seq experiment summaries

Following on from above, **mtab1625** now contains a *SimpleList* object with a single entry named **rnaseq**. For RNA-seq experiments, this entry is a *RangedSummarizedExperiment* object (see package [SummarizedExperiment](#)).

```
sumexp <- mtab1625$rnaseq  
sumexp
```

```
## class: RangedSummarizedExperiment  
## dim: 91080 18  
## metadata(4): pipeline filtering mapping quantification  
## assays(1): counts  
## rownames(91080): EP10SAG000000000001 EP10SAG000000000002 ...  
## OS12G0641500 OS12G0641600  
## rowRanges metadata column names(0):  
## colnames(18): ERR266221 ERR266222 ... ERR266237 ERR266238  
## colData names(9): AtlasAssayGroup organism ... growth_condition  
## sampling_time
```

The matrix of raw counts for this experiment is stored in the *assays* slot:

```
head( assays( sumexp )$counts )
```

```
##                ERR266221 ERR266222 ERR266223 ERR266224 ERR266225  
## EP10SAG000000000001      0         0         0         0         1  
## EP10SAG000000000002      0         0         0         0         0  
## EP10SAG000000000003      0         0         0         0         0  
## EP10SAG000000000004      0         0         0         0         0  
## EP10SAG000000000005      0         0         0         0         0  
## EP10SAG000000000006      0         0         0         0         0  
##                ERR266226 ERR266227 ERR266228 ERR266229 ERR266230  
## EP10SAG000000000001      0         0         0         0         1  
## EP10SAG000000000002      0         0         0         0         0  
## EP10SAG000000000003      0         0         0         0         0  
## EP10SAG000000000004      0         0         0         0         0
```

```

## EP10SAG000000000005      0      0      0      0      0
## EP10SAG000000000006      0      0      0      0      0
## ERR266231 ERR266232 ERR266233 ERR266234 ERR266235
## EP10SAG000000000001      0      1      0      1      1
## EP10SAG000000000002      0      0      0      0      0
## EP10SAG000000000003      0      0      0      0      0
## EP10SAG000000000004      0      0      0      0      0
## EP10SAG000000000005      0      0      0      0      0
## EP10SAG000000000006      0      0      0      0      0
## ERR266236 ERR266237 ERR266238
## EP10SAG000000000001      0      0      0
## EP10SAG000000000002      0      0      0
## EP10SAG000000000003      0      0      0
## EP10SAG000000000004      0      0      0
## EP10SAG000000000005      0      0      0
## EP10SAG000000000006      0      0      0

```

The sample annotations can be found in the *colData* slot:

```
colData( sumexp )
```

```

## DataFrame with 18 rows and 9 columns
##      AtlasAssayGroup      organism      cultivar
##      <character>      <character> <character>
## ERR266221      g5 Oryza sativa Japonica Group Nipponbare
## ERR266222      g2 Oryza sativa Japonica Group Nipponbare
## ERR266223      g2 Oryza sativa Japonica Group Nipponbare
## ERR266224      g5 Oryza sativa Japonica Group Nipponbare
## ERR266225      g3 Oryza sativa Japonica Group Nipponbare
## ...      ...      ...      ...
## ERR266234      g3 Oryza sativa Japonica Group Nipponbare
## ERR266235      g4 Oryza sativa Japonica Group Nipponbare
## ERR266236      g4 Oryza sativa Japonica Group Nipponbare
## ERR266237      g4 Oryza sativa Japonica Group Nipponbare
## ERR266238      g6 Oryza sativa Japonica Group Nipponbare
##      developmental_stage      age
##      <character> <character>
## ERR266221 seedling, two leaves visible, three leaves visible      2
## ERR266222 seedling, two leaves visible, three leaves visible      2
## ERR266223 seedling, two leaves visible, three leaves visible      2
## ERR266224 seedling, two leaves visible, three leaves visible      2
## ERR266225 seedling, two leaves visible, three leaves visible      2
## ...      ...      ...
## ERR266234 seedling, two leaves visible, three leaves visible      2
## ERR266235 seedling, two leaves visible, three leaves visible      2
## ERR266236 seedling, two leaves visible, three leaves visible      2
## ERR266237 seedling, two leaves visible, three leaves visible      2
## ERR266238 seedling, two leaves visible, three leaves visible      2
##      time_unit      organism_part
##      <character>      <character>
## ERR266221      week shoot axis, vascular leaf
## ERR266222      week shoot axis, vascular leaf
## ERR266223      week shoot axis, vascular leaf

```

```
## ERR266224      week shoot axis, vascular leaf
## ERR266225      week shoot axis, vascular leaf
## ...           ...
## ERR266234      week shoot axis, vascular leaf
## ERR266235      week shoot axis, vascular leaf
## ERR266236      week shoot axis, vascular leaf
## ERR266237      week shoot axis, vascular leaf
## ERR266238      week shoot axis, vascular leaf
##              growth_condition sampling_time
##              <character>    <character>
## ERR266221 300 millimolar sodium chloride      5
## ERR266222      normal watering                5
## ERR266223      normal watering                5
## ERR266224 300 millimolar sodium chloride      5
## ERR266225      normal watering                24
## ...           ...
## ERR266234      normal watering                24
## ERR266235 300 millimolar sodium chloride      1
## ERR266236 300 millimolar sodium chloride      1
## ERR266237 300 millimolar sodium chloride      1
## ERR266238 300 millimolar sodium chloride      24
```

Information describing how the raw data files were processed to obtain the raw counts matrix are found in the *metadata* slot:

```
metadata( sumexp )
```

```
## $pipeline
## [1] "iRAP version 0.7.0p1 (http://nunofonseca.github.io/irap/)"
##
## $filtering
## [1] "Discard reads below minimum quality threshold"
## [2] "Check of bacterial contamination; discard offending reads"
## [3] "Discard reads with common uncalled characters (e.g. N)"
## [4] "Remove reads from pair-end libraries that were orphaned by filtering steps 1-3"
##
## $mapping
## [1] "Against genome reference (Ensembl Plants release: 26) tophat2 version: 2.0.12"
##
## $quantification
## [1] "htseq2 version: 0.6.1p1"
```

## Single-channel microarray experiments

Data from a single-channel microarray experiment, e.g. [E-MTAB-1624](#), is represented as one or more *ExpressionSet* object(s) in the SimpleList that is downloaded. *ExpressionSet* objects are indexed by the ArrayExpress accession(s) of the microarray design(s) used in the original experiment.

```
names( mtab1624 )
```

```
## [1] "A-AFFY-126"
```

```
affy126data <- mtab1624[[ "A-AFFY-126" ]]  
affy126data
```

```
## ExpressionSet (storageMode: lockedEnvironment)  
## assayData: 57381 features, 18 samples  
##   element names: exprs  
## protocolData: none  
## phenoData  
##   sampleNames: nippon_control_1hr_rep1 nippon_control_1hr_rep2 ...  
##     nippon_salt_5hr_rep3 (18 total)  
##   varLabels: AtlasAssayGroup organism ... sampling_time (9 total)  
##   varMetadata: labelDescription  
## featureData  
##   featureNames: AFFX-BioB-3_at AFFX-BioB-5_at ...  
##     RPTR-0s-XXU09476-1_at (57381 total)  
##   fvarLabels: probeSets  
##   fvarMetadata: labelDescription  
## experimentData: use 'experimentData(object)'  
## Annotation:
```

The matrix of normalized intensity values is in the *assayData* slot:

```
head( exprs( affy126data ) )
```

```
##               nippon_control_1hr_rep1 nippon_control_1hr_rep2  
## AFFX-BioB-3_at           7.869421           8.365278  
## AFFX-BioB-5_at           7.702652           8.020915  
## AFFX-BioB-M_at           7.652985           8.156907  
## AFFX-BioC-3_at           9.219287           9.556873  
## AFFX-BioC-5_at           8.881807           9.186310  
## AFFX-BioDn-3_at          11.617725          11.877930  
##               nippon_control_1hr_rep3 nippon_control_24hr_rep1  
## AFFX-BioB-3_at           8.637034           8.403105  
## AFFX-BioB-5_at           8.407186           8.150995  
## AFFX-BioB-M_at           8.348117           7.998242  
## AFFX-BioC-3_at           9.895652           9.598754  
## AFFX-BioC-5_at           9.553833           9.354739  
## AFFX-BioDn-3_at          12.193324          11.861900  
##               nippon_control_24hr_rep2 nippon_control_24hr_rep3  
## AFFX-BioB-3_at           8.678257           8.456243  
## AFFX-BioB-5_at           8.413489           8.227663  
## AFFX-BioB-M_at           8.358092           8.111307  
## AFFX-BioC-3_at           9.872856           9.699042  
## AFFX-BioC-5_at           9.594959           9.383014  
## AFFX-BioDn-3_at          12.077460          11.959052  
##               nippon_control_5hr_rep1 nippon_control_5hr_rep2  
## AFFX-BioB-3_at           8.348849           8.643520  
## AFFX-BioB-5_at           8.129436           8.374279  
## AFFX-BioB-M_at           7.978514           8.285401  
## AFFX-BioC-3_at           9.588422           9.828320  
## AFFX-BioC-5_at           9.210903           9.512925  
## AFFX-BioDn-3_at          11.853478          12.043559
```

```
##               nippon_control_5hr_rep3 nippon_salt_1hr_rep1
## AFFX-BioB-3_at           8.401530           8.331911
## AFFX-BioB-5_at           8.193307           8.015213
## AFFX-BioB-M_at           8.046037           7.944433
## AFFX-BioC-3_at           9.685030           9.509499
## AFFX-BioC-5_at           9.379879           9.194149
## AFFX-BioDn-3_at          11.952693           11.800154
##               nippon_salt_1hr_rep2 nippon_salt_1hr_rep3
## AFFX-BioB-3_at           8.463545           8.901247
## AFFX-BioB-5_at           8.290420           8.533720
## AFFX-BioB-M_at           8.139875           8.462567
## AFFX-BioC-3_at           9.676649           9.950136
## AFFX-BioC-5_at           9.343052           9.714590
## AFFX-BioDn-3_at          12.043509           12.263983
##               nippon_salt_24hr_rep1 nippon_salt_24hr_rep2
## AFFX-BioB-3_at           8.363449           8.185702
## AFFX-BioB-5_at           8.112024           7.828596
## AFFX-BioB-M_at           8.115197           7.775437
## AFFX-BioC-3_at           9.549797           9.300713
## AFFX-BioC-5_at           9.276916           8.916605
## AFFX-BioDn-3_at          11.866490           11.647577
##               nippon_salt_24hr_rep3 nippon_salt_5hr_rep1
## AFFX-BioB-3_at           8.586542           8.562492
## AFFX-BioB-5_at           8.207703           8.307241
## AFFX-BioB-M_at           8.231843           8.078984
## AFFX-BioC-3_at           9.624600           9.728441
## AFFX-BioC-5_at           9.371714           9.430824
## AFFX-BioDn-3_at          11.908544           12.029436
##               nippon_salt_5hr_rep2 nippon_salt_5hr_rep3
## AFFX-BioB-3_at           8.366250           8.499076
## AFFX-BioB-5_at           8.275820           8.155382
## AFFX-BioB-M_at           8.114156           8.115064
## AFFX-BioC-3_at           9.632023           9.531450
## AFFX-BioC-5_at           9.358933           9.260982
## AFFX-BioDn-3_at          11.832181           11.987857
```

The sample annotations are in the *phenoData* slot:

```
pData( affy126data )
```

```
##               AtlasAssayGroup               organism
## nippon_control_1hr_rep1          g1 Oryza sativa Japonica Group
## nippon_control_1hr_rep2          g1 Oryza sativa Japonica Group
## nippon_control_1hr_rep3          g1 Oryza sativa Japonica Group
## nippon_control_24hr_rep1         g3 Oryza sativa Japonica Group
## nippon_control_24hr_rep2         g3 Oryza sativa Japonica Group
## nippon_control_24hr_rep3         g3 Oryza sativa Japonica Group
## nippon_control_5hr_rep1          g2 Oryza sativa Japonica Group
## nippon_control_5hr_rep2          g2 Oryza sativa Japonica Group
## nippon_control_5hr_rep3          g2 Oryza sativa Japonica Group
## nippon_salt_1hr_rep1             g4 Oryza sativa Japonica Group
## nippon_salt_1hr_rep2             g4 Oryza sativa Japonica Group
## nippon_salt_1hr_rep3             g4 Oryza sativa Japonica Group
```

|                             |  |
|-----------------------------|--|
| ## nippon_salt_24hr_rep1    | g6 Oryza sativa Japonica Group                     |
| ## nippon_salt_24hr_rep2    | g6 Oryza sativa Japonica Group                     |
| ## nippon_salt_24hr_rep3    | g6 Oryza sativa Japonica Group                     |
| ## nippon_salt_5hr_rep1     | g5 Oryza sativa Japonica Group                     |
| ## nippon_salt_5hr_rep2     | g5 Oryza sativa Japonica Group                     |
| ## nippon_salt_5hr_rep3     | g5 Oryza sativa Japonica Group                     |
| ##                          | cultivar   |
| ## nippon_control_1hr_rep1  | Nipponbare   |
| ## nippon_control_1hr_rep2  | Nipponbare   |
| ## nippon_control_1hr_rep3  | Nipponbare   |
| ## nippon_control_24hr_rep1 | Nipponbare   |
| ## nippon_control_24hr_rep2 | Nipponbare   |
| ## nippon_control_24hr_rep3 | Nipponbare   |
| ## nippon_control_5hr_rep1  | Nipponbare   |
| ## nippon_control_5hr_rep2  | Nipponbare   |
| ## nippon_control_5hr_rep3  | Nipponbare   |
| ## nippon_salt_1hr_rep1     | Nipponbare   |
| ## nippon_salt_1hr_rep2     | Nipponbare   |
| ## nippon_salt_1hr_rep3     | Nipponbare   |
| ## nippon_salt_24hr_rep1    | Nipponbare   |
| ## nippon_salt_24hr_rep2    | Nipponbare   |
| ## nippon_salt_24hr_rep3    | Nipponbare   |
| ## nippon_salt_5hr_rep1     | Nipponbare   |
| ## nippon_salt_5hr_rep2     | Nipponbare   |
| ## nippon_salt_5hr_rep3     | Nipponbare   |
| ##                          | developmental_stage                                |
| ## nippon_control_1hr_rep1  | seedling, two leaves visible, three leaves visible |
| ## nippon_control_1hr_rep2  | seedling, two leaves visible, three leaves visible |
| ## nippon_control_1hr_rep3  | seedling, two leaves visible, three leaves visible |
| ## nippon_control_24hr_rep1 | seedling, two leaves visible, three leaves visible |
| ## nippon_control_24hr_rep2 | seedling, two leaves visible, three leaves visible |
| ## nippon_control_24hr_rep3 | seedling, two leaves visible, three leaves visible |
| ## nippon_control_5hr_rep1  | seedling, two leaves visible, three leaves visible |
| ## nippon_control_5hr_rep2  | seedling, two leaves visible, three leaves visible |
| ## nippon_control_5hr_rep3  | seedling, two leaves visible, three leaves visible |
| ## nippon_salt_1hr_rep1     | seedling, two leaves visible, three leaves visible |
| ## nippon_salt_1hr_rep2     | seedling, two leaves visible, three leaves visible |
| ## nippon_salt_1hr_rep3     | seedling, two leaves visible, three leaves visible |
| ## nippon_salt_24hr_rep1    | seedling, two leaves visible, three leaves visible |
| ## nippon_salt_24hr_rep2    | seedling, two leaves visible, three leaves visible |
| ## nippon_salt_24hr_rep3    | seedling, two leaves visible, three leaves visible |
| ## nippon_salt_5hr_rep1     | seedling, two leaves visible, three leaves visible |
| ## nippon_salt_5hr_rep2     | seedling, two leaves visible, three leaves visible |
| ## nippon_salt_5hr_rep3     | seedling, two leaves visible, three leaves visible |
| ##                          | age time_unit organism_part                        |
| ## nippon_control_1hr_rep1  | 2 week shoot axis, vascular leaf                   |
| ## nippon_control_1hr_rep2  | 2 week shoot axis, vascular leaf                   |
| ## nippon_control_1hr_rep3  | 2 week shoot axis, vascular leaf                   |
| ## nippon_control_24hr_rep1 | 2 week shoot axis, vascular leaf                   |
| ## nippon_control_24hr_rep2 | 2 week shoot axis, vascular leaf                   |
| ## nippon_control_24hr_rep3 | 2 week shoot axis, vascular leaf                   |
| ## nippon_control_5hr_rep1  | 2 week shoot axis, vascular leaf                   |
| ## nippon_control_5hr_rep2  | 2 week shoot axis, vascular leaf                   |
| ## nippon_control_5hr_rep3  | 2 week shoot axis, vascular leaf                   |



```
## nippon_salt_1hr_rep1      2      week shoot axis, vascular leaf
## nippon_salt_1hr_rep2      2      week shoot axis, vascular leaf
## nippon_salt_1hr_rep3      2      week shoot axis, vascular leaf
## nippon_salt_24hr_rep1     2      week shoot axis, vascular leaf
## nippon_salt_24hr_rep2     2      week shoot axis, vascular leaf
## nippon_salt_24hr_rep3     2      week shoot axis, vascular leaf
## nippon_salt_5hr_rep1      2      week shoot axis, vascular leaf
## nippon_salt_5hr_rep2      2      week shoot axis, vascular leaf
## nippon_salt_5hr_rep3      2      week shoot axis, vascular leaf
##                               growth_condition sampling_time
## nippon_control_1hr_rep1      normal watering           1
## nippon_control_1hr_rep2      normal watering           1
## nippon_control_1hr_rep3      normal watering           1
## nippon_control_24hr_rep1     normal watering          24
## nippon_control_24hr_rep2     normal watering          24
## nippon_control_24hr_rep3     normal watering          24
## nippon_control_5hr_rep1      normal watering           5
## nippon_control_5hr_rep2      normal watering           5
## nippon_control_5hr_rep3      normal watering           5
## nippon_salt_1hr_rep1      300 millimolar sodium chloride      1
## nippon_salt_1hr_rep2      300 millimolar sodium chloride      1
## nippon_salt_1hr_rep3      300 millimolar sodium chloride      1
## nippon_salt_24hr_rep1     300 millimolar sodium chloride      24
## nippon_salt_24hr_rep2     300 millimolar sodium chloride      24
## nippon_salt_24hr_rep3     300 millimolar sodium chloride      24
## nippon_salt_5hr_rep1      300 millimolar sodium chloride      5
## nippon_salt_5hr_rep2      300 millimolar sodium chloride      5
## nippon_salt_5hr_rep3      300 millimolar sodium chloride      5
```

A brief outline of how the raw data was normalized is in the *experimentData* slot:

```
preproc( experimentData( affy126data ) )
```

```
## $normalization
## [1] "RMA using oligo (http://www.bioconductor.org/packages/release/bioc/html/oligo.html) version 1.2"
```

## Downloading a single Expression Atlas experiment summary

You can also download data for a single Expression Atlas experiment using the `getAtlasExperiment()` function:

```
mtab3007 <- getAtlasExperiment( "E-MTAB-3007" )
```

```
## Downloading Expression Atlas experiment summary from:
## http://wwwdev.ebi.ac.uk/gxa/experiments/E-MTAB-3007/E-MTAB-3007-atlasExperimentSummary.Rdata ...
## Successfully downloaded experiment summary object for E-MTAB-3007
```